

Bunsen burner

Item no. 1000

DIN-DVGW reg. no. NG-2211AN0133

- air regulation



DIN / Approval

Article no.	Gas
1000	Natural gas
1001	Bottle gas

Notes

Burners according to [DIN 30665](#) cannot be switched off completely. A small flame remains burning even if needle valve or stop cock are turned off. The gas flow shall be disconnected by means of the central stop cock in the laboratory.

Learn more about our laboratory burners [here](#).

The picture shows the burner with safety base (extra charge).

Technical data

Item no.	Max. Temp. (°C)	Consumption (l/h)	Height (mm)	Weight (g)
1000	1500	110	155	210
1001	1540	55	155	210

Item no.	Pipe ø	Head Ø	Bel. (kg)	DIN	GS	Natural gas	Propane	Safety base
1000	14	17	1,25	+	+	+	-	o
1001	14	17	1,35	+	-	-	+	o