Rubber belt clamp



Item no. 2825

The strap clamp is used for the secure fixation of beakers, flasks, bottles and other vessels on the laboratory bench. For this purpose, the stem of the strap clamp is attached, for example, to a laboratory support rod using a bosshead. The vessel to be secured remains standing on the bench and is held in place by the adjustable strap to prevent slipping or tipping over.

The strap clamp is particularly useful when stirrers, shakers, vibrating equipment or connected tubing could cause a vessel to move from its position. It keeps the vessel reliably in place.

- cam type closure
- cork lined jaws
- capacity Ø 60 –170 mm
- rod Ø 10 mm

Technical data

Article no.	Weight
2825	280 g